dental diamant dental diamond

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Diamonds made by Fis - the alternatives

We have been manufacturing diamond tools for dental applications as well as diamond and CBN tools for industrial use since more than 80 years.

Shapes, sizes and abrasive layers of our products are very carefully matched for optimum use in all applications. Custom products have been manufactured for many of our customers. This has resulted in "smart" solutions which are unparalleled in solving technical requirements, offering high quality products.

»High-tech is a standard«

- For dental applications, our production ranges from FG burs for dental offices to plaster trimmer wheels for lab use.
- We manufacture our FG diamonds in varying grit sizes between 504 and 554 (acc. to ISO standards), with choice of one of two types of grit: natural or synthetic (S) diamond (limited to certain sizes and types). The natural diamond crystals typically have an irregular shape and sharp cutting edges which ensure a quick material removal. Layers containing the blocktype, synthetic diamond crystals provide a professional preparation with less surface roughness. Those FG diamonds have a black coating to increase the optical contrast between layer and tooth.
- Only hardened, stainless chromium/nickel steels are used as basic material for the "fis" diamond tools (except a few cutting discs for the lab).
- The tolerances of shaft diameters not only meet the requirements of the ISO standard, but exceed them. With quality in mind, we grind the shafts to close tolerances for precise fit in the handpiece chuck which minimizes runout error and protects the chuck.
- In addition, we have introduced a new generation of diamonds (HighQPlus). These multilayer FG tools provide long life in rough cutting applications.

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General information when using rotating dental diamond grinding points (tools with FG and right angle type shank)

The following instructions provide best results in preparations through careful handling of the diamond tools and their application.

Consideration of the teeth as well as perfect hygienic conditions are the main concern.

- **Die Dental-Diamantschleifstifte müssen vor dem ersten Einsatz sterilisiert werden.**
- **With regard to cleaning/sterilization, the precise, detailed "Directions for use of rotary diamond instruments (FG and right angle type)" in the download area of this website apply.**
- Instruments are to be selected by the head diameter to prevent from damaging neighbour teeth or tissue.
- Rounded shapes should be favoured.
- Burs with coarse grit are not always indicated; usually they require a finish grinding operation.
- Trays or original packing should be used for storing.
- The shafts have to be secured properly and checked for minimum runout before use.
- The burs should be set to the working speed according to the graphical diagram. It should be set before applying the bur to the object.
- Spray cooling with a minimum flow of 50 ml/min must be provided.
- Additional cooling is required when burs longer than 21 mm or heads bigger than 2 mm are being used.
- Intermittent („dabbing“) work allows improved cooling of the grinding zone
- Forces should not exceed 2 N (appr. 200 p). Higher forces do not result in a quicker material removal, but increase the thermal load on the tooth.
- Worn or damaged burs must be replaced immediately.
- Burs should be cleaned in a „cleaning solution“ first after use. See manufacturer's instructions for details.
- Burs can be brushed to remove material sticking to the diamond layer.
- Dental diamond instruments are made of stainless materials. However, they must not get in contact with acids.
- Our dental diamond tools are manufactured from high-quality components providing a long life time. However, we recommend not to use the instruments after an expiration time of 10 years.

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

Cutting discs with penetrated rim

In addition to proven cutting discs with surface layer for all grinding and cutting operations, we are the only manufacturer of cutting discs with penetrated rim covering the diameter range from 10 to 200 mm. This layer is - similar to a sintered bond - completely penetrated with diamond crystals. The discs allow thinnest cuts and extremely long service. They are constantly flexible and appropriate for technicians during cutting and lateral grinding of ceramics, porcelain, plaster and metals.

Directions for use of rotary diamond discs with penetrated rim (Lab tools with electroplated bond, dispersed, Fig. 354).


Cutting discs with electroplated bond and penetrated rim provide unique properties. Other than discs with a single layer, this type achieves a very long life time by a specific multilayer design. Though they are thin and flexible, much more than the ones with sintered bond. They surpass all other types of discs when handled properly.

- The tools should be set to the working speed according to the graphical diagram.
- Cutting discs with penetrated rim are designed for separating. An oscillating motion should be applied while moving in radial direction.
- Although these discs are flexible, lateral grinding should be performed at very low forces only. High lateral forces may result in breakage, as the thickness of the blade usually is only 0,12 mm.
- Bending or turning when cutting into a slot must be avoided.
- A soft cleaning stone is recommended for the cleaning of tools with sintered bond. This procedure will remove the top layer of the bond and release unused diamond grit.

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
Naturdiamantkristalle
Natural diamond crystals

Natural diamond - synthetic diamond

We manufacture our FG tools in varying grit sizes between 504 and 554 (acc. to ISO standards), with choice of one of two types of grit: natural and/or synthetic ("S") diamond (limited to ISO 524 and certain types). The natural diamond crystals typically have an irregular shape and sharp cutting edges which ensure a quick material removal. Layers containing the block-type, synthetic diamond crystals provide a professional preparation with less surface roughness. Those FG diamonds have a black coating to increase the optical contrast between tool and tooth.

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Finzler, Schrock & Kimmel GmbH